Turon

MedTech

en: Instruction for use	lt: Naudojimo instrukcija
bg: Инструкция за употреба	no: Bruksanvisning
cs: Návod k použití	lv: Lietošanas instrukcija
da: Brugsanvisning	nl: Instructies voor gebruik
de: Gebrauchsanweisung	pl: Instrukcja użytkowania
el: Οδηγίες χρήσης	pt: Instruções de utilização
es: Instrucciones de uso	ro: Instrucțiuni de utilizare
et: Kasutusjuhend	sk: Návod na použitie
fi: Käyttöohjeet	sl: Navodila za uporabo
fr: Instructions d'utilisation	sv: Bruksanvisning
hr: Upute za upotrebu	ru: Инструкция по применению
hu: Használati utasítás	tr: Kullanım talimatı
it: Istruzioni per l'uso	uk: Інструкція із застосування

Turon MedTech AB
Kardanvägen 21 • 43232 SE-Varberg • Sweden
www.turonmedtech.com

SteelBushingsInstructionForUse-01

English	04-05
Български език	06-07
Česky	08-09
Dansk	10-11
Deutsch	12-13
Ελληνικά	14-15
Español	16-17
Eesti	18-19
Suomi	20-21
Français	22-23
Hrvatski	24-25
Magyar	26-27
Italiano	28-29
Lietuvių	32-33
Latviešu	34-35
Nederlands	36-37
Norsk	38-39
Polski	40-41
Português	42-43
Română	44-45
Русский	46-47
Slovenčina	48-49
Slovenščina	50-51
Svenska	52-53
Türkçe	54-55
Українська	56-57

English - Instruction for use Needle guide bushings Stainless steel

Intended use

These bushings are designed to guide a biopsy needle and create a sterile barrier when performing biopsy in mammography equipment.

Instructions for use

- Consult Biopsy equipment manufacturer or Biopsy equipment instructions for use to verify that bushings are compatible with biopsy equipment.
- Bushings shall be fixated in a needle guide.
- Bushings shall be sterilized before use.
- Single use bushings are also available.
www.turonmedtech.com
- If any serious incident occurs in relation to the device, the incident shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

REF	Gauge
040501	14 G
040502	16 G
040503	18 G
040504	20 G
040505	20 G open
040508	13 G
040509	12 G
040510	22 G
040511	18 G open

Example of bushings fixated in needle guide.

Cleaning, disinfection and sterilization instruction for needle guide bushings

Caution

New bushings must be cleaned, disinfected and sterilized prior to first use. Note that there are always risks associated with reprocessing of medical devices. Remember that the instructions to follow only is input to the cleaning and sterilization procedures. You are responsible to validate your processes with your local prerequisites such as methods, chemicals, water quality, autoclave, load pattern etc. before use.

Cleaning and Disinfection

You should always try to clean the bushing as soon as possible after usage to ease the cleaning operation. All blood, body fluids and tissue must be completely removed from the bushing prior to sterilization. If brush is used for cleaning only use nylon or similar, never use steel wool. Follow standard ISO 15883 and national guidelines – validated machine cleaning and disinfection is strongly recommended. Pay extra attention to ensure that cleaning and disinfection fluids as well as rinsing water flows properly through the centre of the bushing since this area could be hard to reach and at the same time is critical considering the bushing function. Avoid saline solution and high chloride concentration due to risk of pitting, stress- and spot corrosion. Use cleaning- and disinfection agents known to be compatible with stainless steel alloy EN 1.4305 Always follow the manufacturer's instruction regarding concentration, temperature and exposure time. Do not use detergents with pH > 10,3. The cleaning and disinfection operation, whether it is manual or machined should be finalized by carefully rinse the bushing with cold or preferably demineralized water. After rinsing it should immediately be dried to prevent residual moisture from causing corrosion.

Packaging

If the bushing is supposed to be stored and/or transported after sterilization it must be wrapped. ISO 17665 compliant sterilization pouch should be used together with a chemical indicator to verify that sterilization parameters were met inside.

Sterilization

Preferred sterilization method is steam sterilization with pre-vacuum. Any sterilization process in the range 121°C for 15 minutes to 134°C for 3 minutes in accordance with ISO 17665 works. Make sure your water quality complies with tolerances specified in EN285 or local equivalents to avoid subsequent corrosion or deposits. Also, carefully follow the instructions given by the sterilizer manufacturer. Use suitable pouch divider when placed in autoclave to assist sterilization and drying. Also make sure eventual heavy wraps is placed in the bottom and lighter pouches on the top shelf in the autoclave to enable moisture to drain out without wetting other items in the load.

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